

1 Setup

For these measurements, we used two different setups, the main difference between the two being the absence, or not, of a lead shielding to perform the measurements. The setup with the shielding is given in Figure 1.




Figure 1: Picture of the original setup (with the lead shielding).

The germanium detector is placed in the bottom of the setup, and the radioactive source is placed above (on the blue table on the picture). The beam is collimated through two bricks of lead containing a hole of 10mm diameter and the source is around 202mm above the top of the cryostat (thus 207mm from the top of the crystal). On the sides of the germanium detector are placed two scintillators that will be used to measure coincidences. The whole setup is shielded by lead bricks in order to avoid noise. In the second setup that we used, the source is held by a metallic arm, roughly at the same distance from the top of the crystal and we removed all the lead bricks to measure the peak-to-Compton ratio and the efficiency of the detector.

For these first measurements, we used two sources, one source of Cobalt 60 to perform the calibration, and a source of Cesium 137. The activities of the sources are, respectively, 370kBq and 16.5kBq. Both setups were used with both sources and the measurements last until the net area of the biggest peak in the spectrum was 100000 counts at least, to get proper statistics.

2 Electronics

3 Analysis

The very first thing to perform is the calibration of the detector. Indeed, the MCA only store values of received counts into different channels (we used 8192 channels for these measurements) which has no physical meaning. Thus, we need to calibrate the detector so that we obtain a spectrum representing the number of count received as a function of the energy. In order to do so, we use a source of Cobalt

60 because it contains two very high peaks resulting from its β^- decay (see its decay scheme Figure 3). We used the first setup for this measurement and associated the channel corresponding to each peak to its energy value. The result is given in the Figure 3. Then we used the equation we found to perform the analysis thanks to a python program (the whole program is given in the Appendix).




Figure 2: Decay scheme of Cobalt 60 into Nickel 60. On the final spectrum, we should see two peaks corresponding to the released photons.




Figure 3: Calibration plot representing energy as a function of channel number